

New firmware version 01.04.14 for your CoaguChek INRange meter

Date:

June 2018

Rev. 1.0

We herewith introduce the new firmware version 01.04.14 for your CoaguChek INRange meter.

Please note:

The difference between the former firmware version 01.04.12 and the new firmware version 01.04.14 is explained as follows:

Due to resolution problems, the Airwave Encryption Key 2d barcode, displayed in the "About" menu, could not be read out with the old firmware version 01.04.12.

The new firmware version 01.04.14 addresses this issue and has a higher resolution now. The 2d barcode can now be read.

Should you have any question about the CoaguChek INRange firmware version 01.04.14, please contact your Roche Service organization.

Sincerely